## NICO BrainPath® Instructions for Use

### Table of Contents

|                  |    |
|------------------|----|
| English .....    | 2  |
| Deutsch .....    | 10 |
| Español .....    | 17 |
| Français .....   | 24 |
| Italiano .....   | 31 |
| Nederlands ..... | 38 |
| Português .....  | 45 |

# ENGLISH

## **NICO BrainPath® Instructions for Use**

**Caution:** Federal (US) Law restricts this device to sale by or on the order of a physician.

Please read all information carefully. Failure to properly follow these instructions may lead to unintended surgical consequences. Important: This package insert is designed to provide instructions for use for the NICO BrainPath. It is not a reference to surgical techniques.

### **Intended Use**

To provide for access and allow for visualization of the surgical field during brain and spinal surgery. Indications may include subcortical access to diseases such as the following:

- Primary/Secondary brain tumors
- Vascular abnormalities/malformations
- Intraventricular tumors/cysts

### **Contraindications**

None known

### **Warnings**

- Use of the NICO BrainPath should be performed only by surgeons having adequate training and experience with using a port to access the brain or spine.
- Avoid use and/or take special care in areas of the brain which carry increased risk.
- Do not resterilize or reuse the NICO BrainPath Sheath or sterile accessories. These sterile devices are intended for a single use only. Resterilization may compromise their integrity, which may result in it malfunctioning.
- The sheath and obturator of the NICO BrainPath which come into contact with bodily fluids require special disposal and handling to prevent biological contamination.
- If prion contamination is known or suspected, it is recommended the reusable components (obturator, manipulation tool and navigation stylet) be incinerated.
- The BrainPath has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of BrainPath in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

### **Precautions**

Devices and accessories manufactured or distributed by companies not authorized by NICO Corporation may not be compatible with the NICO BrainPath. Use of such products may lead to unanticipated results and possible injury to the user or patient.

### **Special Facilities, Training, or Qualifications for Users**

BrainPath is prescription-use-only and intended to be used by a licensed neurosurgeon. There are no special facility, training, or user qualifications.

**Symbols (see primary product label for applicable symbols)**

|                                                                                  |                                                                 |                                                                                   |                                                                             |
|----------------------------------------------------------------------------------|-----------------------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------|
|  | Prescription Use Only                                           |  | Catalog Number                                                              |
|  | Single use only and should not be reused                        |  | Serial Number                                                               |
|  | Caution                                                         |  | Non-Sterile                                                                 |
|  | Indicates that the device is not made with natural rubber latex |  | Indicates that the device has been sterilized via gamma radiation processes |
|  | Indicates the device is non-pyrogenic                           |  | Indicates use-by date                                                       |
|  | Indicates Lot Number                                            |  | Indicates the device manufacturer                                           |
|  | Indicates the manufacture date                                  |  | Do not use if package is damaged                                            |
|  | Consult the IFU                                                 |  | Indicates item is a Medical Device                                          |
|  | Indicates a double sterile barrier                              |  | Unique Device Identifier                                                    |

**Device Description**

The BrainPath consists of multiple-sized reusable and re-sterilizable obturators with coordinating single patient use disposable sheaths. The obturator and sheath are assembled in the operating room immediately prior to use. After placement, the obturator is removed leaving behind the sheath which provides a 13.5 mm or 11 mm surgical corridor.

The obturators are designed with a channel and set-screw for placement and securement of third-party navigation probes. The channel of the obturator stops before the tapered tip begins thus navigation settings must reflect the obturator tip length as an extension/offset. For navigation probes that cannot be placed in the existing channel, various adapters are available that ensure the tip of the probe is aligned with the head of the obturator, thus navigation settings must reflect the entire obturator length as an extension/offset. The following tables outline the available BrainPath configurations and relevant obturator dimensions for navigation:


**Table 1: BrainPath Configurations**

| Description          | Sheath Inner-Diameter (mm) | Sheath Length (mm) | Obturator Tip Length (mm) |
|----------------------|----------------------------|--------------------|---------------------------|
| 13.5 mm x 50 mm (ST) | 13.5                       | 50                 | 7.5                       |
| 13.5 mm x 50 mm      |                            | 50                 | 15                        |
| 13.5 mm x 60 mm      |                            | 60                 |                           |
| 13.5 mm x 75 mm      |                            | 75                 |                           |
| 13.5 mm x 95 mm      |                            | 95                 |                           |
| 11 mm x 50 mm        | 11                         | 50                 |                           |
| 11 mm x 60 mm        |                            | 60                 |                           |
| 11 mm x 75 mm        |                            | 75                 |                           |

Table 2: Relevant Obturator Dimensions

|                      | Third-Party Nav. Probe:<br>Within Obturator Channel |                                                     | Third-Party Nav. Probe:<br>External using NICO Adapter |                                                       |
|----------------------|-----------------------------------------------------|-----------------------------------------------------|--------------------------------------------------------|-------------------------------------------------------|
| Obturator            | Acceptable<br>Diameter<br>(mm)                      | Obturator Tip-<br>Length/ Extension/<br>Offset (mm) | Acceptable<br>Diameter<br>(mm)                         | Obturator Overall<br>Length/Extension/<br>Offset (mm) |
| 13.5 mm x 50 mm (ST) | < 4.3                                               | 7.5                                                 | 4.3 - 5.2                                              | 73                                                    |
| 13.5 mm x 50 mm      |                                                     | 15                                                  |                                                        | 81                                                    |
| 13.5 mm x 60 mm      |                                                     |                                                     |                                                        | 91                                                    |
| 13.5 mm x 75 mm      |                                                     |                                                     |                                                        | 106                                                   |
| 13.5 mm x 95 mm      |                                                     |                                                     |                                                        | 126                                                   |
| 11 mm x 50 mm        |                                                     |                                                     |                                                        | 92                                                    |
| 11 mm x 60 mm        |                                                     |                                                     |                                                        | 102                                                   |
| 11 mm x 75 mm        |                                                     |                                                     |                                                        | 116                                                   |

Figure 1: Obturator Measurement Locations

The BrainPath accessories include a “manipulation tool,” which is similar to a dental probe and is used for manipulating the position of the sheath after it has been placed. The accessories also include “Shepherd’s Hooks” for attaching to various commercially available retractors, a “navigation stylet” for use with the 95 mm obturator, “Probe Adapters” for use with third-party navigation probes that cannot be placed in the channel of the obturator, and a sterilization tray for sterilization of the reusable components (i.e., obturators, manipulation tools, and stylet).

Table 3: NICO BrainPath Components Supplied

| BrainPath Component          | Sterile or Non-Sterile | Single-Use or Reusable |
|------------------------------|------------------------|------------------------|
| Sheath (all sizes)           | Sterile                | Single-Use             |
| Obturators (all sizes)       | Non-Sterile            | Reusable               |
| Manipulation Tool            | Non-Sterile            | Reusable               |
| Navigation Stylet            | Non-Sterile            | Reusable               |
| Sterilization Tray           | Non-Sterile            | Reusable               |
| Shepherd’s Hooks (all types) | Sterile                | Single-Use             |
| Probe Adapter (all types)    | Sterile                | Single-Use             |

## **Instructions for Use**

1. The surgeon determines and plans the appropriate approach to the area of interest and prepares the patient according to standard operative preparation protocols as required by the facility.
2. Prior to use, inspect packaging of the sheath and/or sterile accessory (e.g., Shepherd's Hook) for the appropriate size as requested by the surgeon and for any damage or gross deformations. If there is any indication that sterility or functionality has been compromised, obtain a new device and contact NICO Customer Service about returning the original.
3. Remove sheath from packaging using sterile technique.
4. Obtain sterile obturator and confirm that size of obturator matches size of sheath as requested by the surgeon.
5. Obtain sterile manipulation tool and/or navigation stylet (if desired).
6. Obtain appropriate sterile Shepherd's Hook (if desired).
7. Engage the Shepherd's Hook with the proximal portion of the sheath (if desired) just below the knurled ring.
8. Wet the shaft of obturator using sterile saline solution.
9. Wet the entire sheath by dipping it in sterile saline solution
10. Properly assemble obturator and sheath. Ensure that sheath is fully seated on obturator such that there is a smooth transition between the distal end of the sheath and obturator.
11. If desired, place the third-party navigation probe within the obturator channel and secure using thumb screw. Alternatively, for navigation probes that cannot be placed in the obturator channel, obtain the appropriate probe adapter, secure, and then place the third-party navigation probe into the adapter. For the 95 mm Obturator, the 95 mm navigation stylet may be used as a post for attachment of a clamping third-party optical navigation component (if desired).
12. After the creation of an appropriate craniotomy and dural incision, insert BrainPath assembly into tissue to desired location. NOTE: Depth bands exist on obturator and sheath and may be used as a reference during placement.
13. Once the desired location has been reached, advance sheath until the most proximal white band on the obturator appears and is immediately adjacent to the knurled ring on the sheath. Once advanced this distance, the distal end of the sheath will be aligned with the tip of the obturator.
14. While holding the sheath in the desired location, remove obturator from the sheath and the surgical field. Once obturator has been removed, sheath position may be manipulated using manipulation tool. Sheath may also be secured further by placing sutures through the holes present in the knurled ring of the sheath.
15. If using a Shepherd's Hook, secure the proximal end of the hook into the applicable retractor arm.
16. If using a Select Lock Shepherd's Hook, secure the rotational position of the sheath (if desired) by tightening the distal tapered nut with the textured surface.

17. Perform procedure using surgeon desired instrumentation through the BrainPath sheath as appropriate.
18. Once the surgeon has completed the desired end result through the BrainPath, remove sheath and dispose per appropriate hospital procedure.

### **How Supplied**

The BrainPath Sheath and accessories are supplied sterile and packaged for single use. Discard into an appropriate container after use. The BrainPath Obturator, Manipulation Tool and Stylet are packaged non-sterile and require onsite processing/sterilization prior to use. Refer to “BrainPath Reprocessing Instructions” (LL-5904) for detailed information concerning onsite reprocessing.

### **Special Storage or Handling Conditions**

There are no special storage or handling conditions for the BrainPath single-use components. For reusable components, refer to the BrainPath Reprocessing instructions.

### **Device Related Serious Incident Reporting**

Any serious incident or serious adverse event involving these products should be reported to NICO or the local Distribution Partner immediately. Serious incidents may also need to be reported to the local Regulatory Authority (e.g., Competent Authority of the Member State in EU). To report a serious incident, contact the local Representative or the Customer Service at 317-660-7118 ext.100.

### **Clinical Benefits of BrainPath**

Published literature on procedures that used BrainPath suggests the following clinical benefits:

- Reduction in parenchymal injury

### **Clinical Risks**

All surgeries carry risks; the following table outlines the adverse events mentioned in the literature inherent to neurosurgical procedures for which BrainPath could be used. These risks are known risks associated with these neurosurgical procedures and underlying diseases. None of the studies stated the events or deaths were attributed to the use of the BrainPath. Post-Market Surveillance data also continues to support low complaints for the BrainPath.

Procedural Risks Table

| Event           | References                      | Event                   | Reference                       |
|-----------------|---------------------------------|-------------------------|---------------------------------|
| Aphasia         | Mansour et al., 2020            | Iatrogenic hypervolemia | Bauer et al., 2017              |
| Arm pain        | Marenco-Hillebrand et al., 2020 | Ischemia                | Mansour et al., 2020            |
| Blood loss      | Marenco-Hillebrand et al., 2020 | Leg weakness            | Marenco-Hillebrand et al., 2020 |
| Brain contusion | Marenco-Hillebrand et al., 2020 | MCA perforator          | Monroy-Sosa et al., 2021        |

|                                   |                                                                                                                                                |                                                |                                 |
|-----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------|---------------------------------|
| <b>Cardiac demand</b>             | Marenco-Hillebrand et al., 2020                                                                                                                | <b>Memory and word recall impairment</b>       | Eichber et al., 2020a           |
| <b>Cerebrospinal fluid</b>        | Marenco-Hillebrand et al., 2020                                                                                                                | <b>Paresthesias</b>                            | Lin et al., 2020                |
| <b>Death</b>                      | Bauer et al., 2017; Echeverry et al., 2020; Griessenauer et al., 2018; Przybylowski et al., 2015; Rutkowski et al., 2019; Eliyas et al., 2016; | <b>Postoperative stroke with unknown cause</b> | Marenco-Hillebrand et al., 2020 |
| <b>Deep vein thrombosis</b>       | Bakhsheshian et al., 2019; Gassie et al., 2019                                                                                                 | <b>Rebleeding</b>                              | Mansour et al., 2020            |
| <b>Early neurological deficit</b> | Eichberg et al., 2019                                                                                                                          | <b>Reoperation</b>                             | Rutkowski et al., 2019          |
| <b>Encephalomyelitis</b>          | Gassie et al., 2019                                                                                                                            | <b>Seizure</b>                                 | Eichber et al., 2020a           |
| <b>Haemorrhage</b>                | Mansour et al., 2020                                                                                                                           | <b>Stroke</b>                                  | Eichberg et al., 2018           |
| <b>Hematoma</b>                   | Chen et al., 2016                                                                                                                              | <b>Surgical wound dehiscence</b>               | Goren et al., 2019              |
| <b>Hemiparesis</b>                | Bakhsheshian et al., 2019; Eliyas et al., 2016                                                                                                 | <b>Transient short</b>                         | Marenco-Hillebrand et al., 2020 |
| <b>Hemiplegia</b>                 | Akbari et al., 2019; Marenco-Hillebrand et al., 2020; Mansour et al., 2020                                                                     | <b>Upper extremity</b>                         | Marenco-Hillebrand et al., 2020 |
| <b>Hemiplegia / hemiparesis</b>   | Marenco-Hillebrand et al., 2020                                                                                                                | <b>Ventriculitis</b>                           | Rutkowski et al., 2019          |
| <b>Homonymous / hemianopsia</b>   | Lin et al., 2020                                                                                                                               | <b>Worsening of preoperative dysphasia</b>     | Marenco-Hillebrand et al., 2020 |
| <b>Hydrocephalus</b>              | Marenco-Hillebrand et al., 2020                                                                                                                |                                                |                                 |

## Summary of Safety & Clinical Performance

A summary of safety and clinical performance can be found on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

## Basic-UDIs

|                     |                |
|---------------------|----------------|
| Obturator           | 081674402BP1DA |
| Sheaths             | 081674402BP2DC |
| 95 mm Stylet        | 081674402BP3DE |
| Shepherd's Hooks    | 081674402BP4DG |
| Probe Adapter       | 081674402BP5DJ |
| Reusable Components | 081674402BP6DL |

## For More Information

For more information on the NICO BrainPath, contact your NICO representative or NICO Customer Service at 317-660-7118 ext.100. Additional information may be found at [www.niconeuro.com](http://www.niconeuro.com), including electronic versions of the labeling: [www.niconeuro.com/labeling](http://www.niconeuro.com/labeling)

Patent: [www.niconeuro.com/patents](http://www.niconeuro.com/patents)

© 2021 NICO Corporation. All Rights Reserved. Printed in the USA.

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© 2021 NICO Corporation. Todos os direitos reservados.