

# NICO APS Product Information

NICO  
Automated  
Preservation  
System™

## A Systems Solution

*Automated tissue collection  
and biological preservation,  
while in the operating suite*



Scan for  
APS setup  
tutorial

## Automated Preservation System

exclusively for use with the Myriad NOVUS®



### MORE TISSUE

*ALL tissue removed with the Myriad  
is collected for post-procedural analysis*

### BETTER TISSUE


*Immediate refrigeration of tissue,  
while limiting exposure to the atmosphere*

### LESS VARIABILITY


*Standardizes intraoperative specimen  
handling and reduces variability  
that may occur with manual collection*

### ENABLES EFFICIENCIES

*Fully automated system providing  
an easy process for regionalizing samples  
by intratumor location*



✓ Specimen Collection



✓ Specimen Infusion




✓ Specimen Preservation

# Automated Preservation System Ordering Information

| Automated Preservation System: Pack Options   | Catalog SKU # | Quantity |
|---|---------------|----------|
| Packs Include:<br>1-Specimen Collector with filter element and labeling tag, 1- 0.76mm Specimen Infusion Valve and 1-Specimen Preserver |               |          |
| Automated Preservation System Pack with Clamshell Filter  | NN-1010       | Each     |
| Automated Preservation System Pack with Scoop Filter  | NN-1011       | Each     |
| A La Carte Options  | Catalog SKU # | Quantity |
| Specimen Collector with Clamshell Filter  | NN-9002       | Box of 5 |
| Specimen Collector with Scoop Filter  | NN-9008       | Box of 5 |
| Filter Element Replacements - Clamshell   | NN-9003       | Box of 5 |
| Filter Element Replacements - Scoop   | NN-9009       | Box of 5 |
| Specimen Infusion Valve with 0.5mm Metering Line  | NN-8016       | Box of 5 |
| Specimen Infusion Valve with 0.76mm Metering Line   | NN-8039       | Box of 5 |

Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution company, or customer service for the latest product offerings.



The NICO products identified herein are covered by one or more of the following U.S. patents: US8357175, US8430825, US8460327, US8496599, US8657841, US8702738, US8888803, US9028518, US9216031, US9445831, US9504247, US5911701, US7556622, US9279751, US6245084, as well as other patents pending worldwide. Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution or company, customer service for the latest product offerings LIT461RevB

