

BrainPath® Product Information

NICO
BrainPath®

*Minimally-disruptive, navigable,
trans-sulcal access through
the subcortical space*



BrainPath supports
Minimally Invasive
Parafascicular Surgery
(MIPS)


BrainPath is part of a systems approach for parafascicular, trans-sulcal access to address:

- Primary and secondary tumors
- Vascular abnormalities and malformations
- Intraventricular tumors and cysts



Product Features

- Atraumatic tip designed to minimize tissue disruption during advancement to the surgical site
- Securely locks in navigation pointers during placement
- Frosted outer sheath designed to visualize nearby anatomy and minimize light reflection
- Variable obturator lengths to reach shallow or deep subcortical abnormalities
- Two obturator diameters, including the world's smallest option for trans-sulcal access, enabling bi-manual surgical technique





Scan the QR Code
to View Clinical and
Economic Evidence
using BrainPath

NICO BrainPath Ordering Information

Capital Components	Catalog SKU #	Quantity
Kit Includes: 2 Obturators of each size/diameter, 4 Manipulation tools, 2 Sterilization trays, 2 BrainPath stylets (95mm)		
BrainPath Family Capital Kit	NN-7044	Each
Access: Pack Options	Catalog SKU #	Quantity
Packs Include: 1 BrainPath Sheath, 1 pre-bent Myriad Handpiece, 1 NICO meniGLIDE, 1 set of SecureHold Retractors, 1 Select-Lock Greenberg Shepherd's Hook		
BrainPath Procedure Pack (11x50 sheath with 13GA Myriad)	NN-2000	Each
BrainPath Procedure Pack (11x60 sheath with 13GA Myriad)	NN-2001	Each
BrainPath Procedure Pack (11x75 sheath with 13GA Myriad)	NN-2002	Each
BrainPath Procedure Pack (13.5x50 sheath with 13GA Myriad)	NN-2003	Each
BrainPath Procedure Pack (13.5x60 sheath with 13GA Myriad)	NN-2004	Each
BrainPath Procedure Pack (13.5x75 sheath with 13GA Myriad)	NN-2005	Each
BrainPath Procedure Pack (13.5x95 sheath with 13GA Myriad)	NN-2006	Each
BrainPath Procedure Pack (13.5x50 sheath with 11GA Myriad)	NN-2024	Each
BrainPath Procedure Pack (13.5x60 sheath with 11GA Myriad)	NN-2025	Each
BrainPath Procedure Pack (13.5x75 sheath with 11GA Myriad)	NN-2026	Each
BrainPath Procedure Pack (13.5x95 sheath with 11GA Myriad)	NN-2027	Each
Access - Illumination - Removal: Pack Options		
Packs Include: 1 BrainPath Sheath, 1 Pre-Bent Myriad Handpiece, 1 NICO meniGLIDE, 1 set of SecureHold Retractors, 1 Select-Lock Greenberg Shepherd's Hook, 1 Myriad-LX Fiber Optic Cable, 2 Myriad-LX Sleeves		
BrainPath + Illumination Pack (11x50 sheath with 13GA Myriad)	NN-3000	Each
BrainPath + Illumination Pack (11x60 sheath with 13GA Myriad)	NN-3001	Each
BrainPath + Illumination Pack (11x75 sheath with 13GA Myriad)	NN-3002	Each
BrainPath + Illumination Pack (13.5x50 sheath with 13GA Myriad)	NN-3003	Each
BrainPath + Illumination Pack (13.5x60 sheath with 13GA Myriad)	NN-3004	Each
BrainPath + Illumination Pack (13.5x75 sheath with 13GA Myriad)	NN-3005	Each
BrainPath + Illumination Pack (13.5x95 sheath with 13GA Myriad)	NN-3006	Each
BrainPath + Illumination Pack (13.5x50 sheath with 11GA Myriad)	NN-3007	Each
BrainPath + Illumination Pack (13.5x60 sheath with 11GA Myriad)	NN-3008	Each
BrainPath + Illumination Pack (13.5x75 sheath with 11GA Myriad)	NN-3009	Each
BrainPath + Illumination Pack (13.5x95 sheath with 11GA Myriad)	NN-3010	Each
A La Carte Options	Catalog SKU #	Quantity
Myriad Canister	NN-9001	Box of 10
Select-Lock Shepherd's Hook - Greenberg	NN-8042	Box of 2
Shepherd's Hook - Greenberg (original)	NN-8024	Box of 2
Shepherd's Hook - Sugita	NN-8025	Box of 2
Shepherd's Hook - Budde	NN-8026	Box of 2
BrainPath Navigation Probe Adapter	NN-8045	Box of 2

Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution company, or customer service for the latest product offerings.



The NICO products identified herein are covered by one or more of the following U.S. patents: US8357175, US8430825, US8460327, US8496599, US8657841, US8702738, US8888803, US9028518, US9216031, US9445831, US9504247, US5911701, US7556622, US9279751, US6245084, as well as other patents pending worldwide. Not all products shown may be available in all regulatory jurisdictions. Consult your local representative, distribution or company, customer service for the latest product offerings. LIT208RevP

NICO
CORPORATION